# bosotherm xpress

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# bosotherm xpress

## Digital infrared-forehead-thermometer

### Instruction manual

The bosotherm xpress infrared-forehead-thermometer is a high quality product incorporating state of the art technology. With its unique technology the bosotherm xpress provides a stable reading with each measurement. The instrument performs a self-test every time it is switched on. This ensures the specified accuracy of measurements.

The bosotherm xpress is intended for the monitoring of human body temperature of persons of all ages.

Please read these instructions carefully before using the instrument and keep them in a safe place.

Table of Contents	Page
<b>1. The Advantages of your bosotherm xpress</b>	<b>17</b>
<b>2. Important Safety Instructions</b>	<b>18</b>
<b>3. Product Description</b>	<b>18</b>
<b>4. How the bosotherm xpress measures forehead temperature</b>	<b>19</b>
<b>5. How to use the bosotherm xpress</b>	<b>19</b>
<b>5.1 Important notes for the valuation of the measurements</b>	<b>20</b>
<b>6. How to recall stored values</b>	<b>21</b>
<b>7. Changing from Fahrenheit to Celsius and vice-versa</b>	<b>22</b>
<b>8. Control Displays and Symbols</b>	<b>23</b>
<b>9. Error Messages</b>	<b>24</b>
<b>10. Cleaning and Storage</b>	<b>25</b>
<b>11. Calibration</b>	<b>26</b>
<b>12. Replacing the Battery</b>	<b>26</b>
<b>13. Warranty</b>	<b>26</b>
<b>14. Care for the environment</b>	<b>27</b>
<b>15. Technical Specifications</b>	<b>28</b>

## **1. The Advantages of your bosotherm xpress**

---

### **Measurement in 1 second**

The innovative infrared technology allows measurement of the forehead temperature in only 1 second.

### **Accurate and reliable**

Due to the unique construction of the probe, the advanced infrared sensor, and the complete calibration process, this unit offers a very accurate and reliable temperature measurement comparable to measurements taken from the oral site with a regular digital thermometer.

### **Comfortable and easy handling**

- No probe covers are required. This makes daily use of the thermometer simple and cost effective.
- Measurements with the bosotherm xpress are possible even while a child is sleeping.

### **Auto-Display Memory**

The last reading is automatically displayed for 3 seconds when the unit is switched on ("M" is displayed simultaneously). Storage and recall of the last 12 readings enabling efficient tracking of temperature variations.

### **Safe and Hygienic**

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.

### **Fever Alarm**

10 short beeps and the flashing of the measured temperature alert the patient that the temperature is higher than 37.5 °C (99.5 °F) and she/he may have fever.

## **2. Important Safety Instructions**

---

- Never use the thermometer for purposes other than those it has been intended for. Please follow the general safety precautions. Details regarding special safety precautions are available from the supplier on demand.
- Never immerse the bosotherm xpress into water or other liquids as it is not waterproof. For cleaning and disinfecting please follow the instructions in the «Cleaning and Storage» section.
- Do not expose the instrument directly to the sun and keep it in a dust-free and dry area.
- Do not use the thermometer if there are signs of damage on the measuring tip or on the instrument itself. If damaged, do not attempt to repair the instrument by yourself! Please contact your nearest authorized service centre or contact boso directly.
- The bosotherm xpress consists of high-quality precision parts. Do not drop the instrument! Protect it from severe impact and shock. Do not twist the instrument and the measuring probe!
- Electrical medical devices have to comply with the special safety precautions according to the electromagnetical toleration. According notes are available from the manufacturer upon request.

## **3. Product Description**

---

**(1) Probe**

**(2) Start button**

**(3) LCD Display**

**(4) O/I button**

**(5) Battery Cover**


**(1)**

**(2)**

**(3)**

**(4)**

**(5)**

#### 4. How the bosotherm xpress measures forehead temperature

The bosotherm xpress measures infrared energy radiated from the artery under the skin of the forehead and the surrounding tissue. This energy is collected through the lens and converted to a temperature value. **The most accurate result is obtained when measuring the temperature directly at the centre of the forehead.**

#### 5. How to use the bosotherm xpress

**IMPORTANT:** Before taking the first measurement, please remove the safety slip from the battery slot.

1. Press the O/I button. The LCD is activated and shows all segments for 2 seconds.
2. The last measurement reading will be shown on the display automatically for 3 seconds together with the «M» icon.
3. When the °C or °F icon is flashing, a beep sound is heard and the thermometer is ready for the measurement.
4. Place the probe at the centre of the forehead. Make sure that the tip of the probe is completely in contact with the skin of the forehead.
5. Press the «START» button briefly. A beep sound confirms the end of the measurement.
6. Remove the thermometer from the forehead. The LCD displays the measured temperature.

**NOTE:** 10 short beeps and the flashing of the measured temperature alert the patient that the temperature is higher than 37.5 °C (99.5 °F) and she/he may have fever.

7. Please wait until the °C- or °F symbol flashes before measuring again. If so, continue like explained from 4. on.
8. The thermometer switches off automatically about 1 minute after the last measurement.

## 5.1 IMPORTANT NOTES FOR THE VALUATION OF THE MEASUREMENTS

---

Since body temperature varies from person to person, it's useful to know the individual range of the own body temperature.





To get to know this individual range, as well as to get to know how the infrared thermometer works, we recommend to do some test measurements. These test measurements should be done when you are in a healthy, fever free condition.

**The standard range of the body temperature for oral and forehead measurements is 35,7 – 37,2 °C, for rectal measurements it is 36,1 – 37,6 °C.**

- Always take the temperature in the same spot, since the temperature readings may vary at different spots.
- The thermometer as well as the patient should stay in the same ambient temperature for at least 10 minutes before the measurement takes place.
- Readings could be too high, if hair covers the measuring spot.
- After sleeping please wait for approximately 5 minutes before taking the forehead temperature.
- Please do not use the thermometer in places with high humidity.
- Before measuring the temperature, please remove sweat or make up from the forehead.
- It is possible that the forehead temperature differs from the body temperature (which for example was measured rectally) due to a contraction of the blood vessels, which can occur in every fever stage. In this case, the skin feels cold and the temperature measured with the infrared forehead thermometer could be noticeably low.
- Skin temperature of babies reacts quickly to the ambient temperature. Please do not measure the temperature while or immediately after changing the nappy; the measured temperature could be too low.

**IN CASE OF DOUBT A COMPARATIVE MEASUREMENT IS RECOMMENDED.**

## 6. How to recall stored value



Display	Display Meaning	Description
	Recall Mode	Press the START button to enter Recall Mode when the unit is switched off. The memory icon « <b>M</b> » flashes.
 	Reading 1 (the last reading)	Press and release the START button to recall the last reading. « <b>1</b> » appears on the display together with the memory icon « <b>M</b> ». Now the last reading appears.
 	Reading 2 (the reading before last)	Press and release the START button again to recall the reading before last.
 	Reading 12 (the oldest reading)	Press and release the START button consecutively to recall readings in succession, up to the last 12 readings.

Pressing and releasing the START button after the last 12 readings will resume the above sequence from reading 1.






## 7. Changing from Fahrenheit to Celsius and vice-versa

---







The bosotherm xpress Forehead Thermometer can display temperature measurements in either Fahrenheit or Celsius. To change the temperature scale from Fahrenheit to Celsius or vice versa, switch off the thermometer and press START for about 10 seconds until a flashing °F or °C icon is shown on the LCD. Press and release START to switch from °F to °C, or vice versa. After choosing the correct scale wait for about 5 seconds until the thermometer changes automatically into the measuring mode.

## 8. Control Displays and Symbols

Display	Display Meaning	Description
	All segments displayed	Press the O/I button to switch on the thermometer, all segments will be shown for 2 seconds.
 	Memory	The last reading will be shown automatically on the display for 3 seconds.
 	Ready	The unit is ready for the measurement when the °C or °F icon is flashing.
 	Measurement completed	The reading appears on the display while the °C or °F icon is shown permanently. The thermometer is ready for measuring again when the °C or °F icon is flashing.
	Battery change	The battery should be changed when the battery symbol appears on the display while the thermometer is switched on.

## 9. Error Messages

Display	Display Meaning	Possible cause and fault correction
 	Measured temperature too high	Thermometer displays «H» when the measured temperature is higher than 42.2°C or 108.0°F.
 	Measured temperature too low	Thermometer displays «L» when the measured temperature is lower than 34.0°C (93.2°F).
 	Ambient temperature too high	Thermometer displays «H» together with the «▲» symbol when the ambient temperature is higher than 40.0 °C or 104.0 °F.
 	Ambient temperature too low	Thermometer displays «L» together with the «▼» symbol when the ambient temperature is lower than 16.0 °C or 60.8 °F.
	Malfunction display (error)	Please contact your local service centre or bosso if the display shows the error sign consecutively.
	Blank display	Please check if the battery is in a good condition. Also check the polarity (<+> and <->) of the battery.
	Dead battery indication	The battery should be replaced immediately if the battery icon is the only symbol shown on the display.

## **10. Cleaning and storage**

Use an alcohol swab or cotton swab moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer.

**Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids.**

Take care not to scratch the surface of the LCD. Remove the battery from the instrument if it is not required for extended periods of time in order to avoid damage to the thermometer resulting from a leaking battery.

**The measuring probe must be completely dry before a new measurement takes place.**

**Correct measurements are only guaranteed if the measuring probe is clean and not damaged in any way.**

## 11. Calibration

---

Regular checks on the accuracy of the measurement system should be carried out in accordance with the valid directives of the individual countries, governing this issue.

## 12. Replacing the Battery

---

The bosotherm xpress is supplied with one lithium battery, type CR2032.

Replace with a new CR2032 battery when the battery symbol appears flashing or permanently on the LCD display.

To open the battery slot, press with the thumb on the marked spot and slide open the lid downwards. Change the battery (check correct polarity) and close the battery slot by sliding the lid upwards.

## 13. Warranty

---

For details on the conditions of warranty, please refer to the card "Certificate of Warranty" supplied. A warranty claim must be submitted within the warranty period. For both guarantee and repair work, please bring or send the defective thermometer carefully packed and postage paid either to your authorized dealer or directly to:

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Enclose either the "Certificate of Warranty" completed and stamped by the dealer or attach the sales slip/invoice.

## 14. Care for the environment

---

Used batteries and defective electronic devices do not belong into the household waste.

They can be deposited at special centres. Please contact your local authorities for information.

## 15. Technical Specifications

<b>Type:</b>	<b>Digital infrared forehead thermometer bosotherm xpress</b>
Measuring Range:	34.0 °C to 42.2 °C (93.2 °F to 108.0 °F)
Accuracy under laboratory conditions:	±0.2 °C, 34.0 ~ 42,2 °C (±0.4 °F, 93,2 ~ 108,0 °F)
Clinical repetition for:	small infants: 0,09 °C Infants: 0,08 °C Adults: 0,1 °C
Display:	Liquid Crystal Display, 4 digits, 0,1 °C scale (0,1°F scale)
Memory:	Recall of 12 readings.
Display Backlight:	a. The display will be lighted for 4 seconds when the unit is turned ON. b. The display will be lighted again for 5 seconds when the measurement has been completed
Operating temperature:	16°C to 40 °C (60.8 °F to 104 °F)
Storage/transport temperature:	-20 °C to +50 °C (-4 °F to 122 °F)
Automatic-Switch-off:	Approx. 1 minute after last measurement has been taken
Battery:	CR2032 3 V BATTERY (X1) - at least 1000 measurements.
Classification:	Type: BF, (  )
Dimensions:	120 mm (L) x 50 mm (W) x 50mm (H)
Weight:	49g (with battery), 45g (w/o battery)
Standards:	Complies with ASTM E-1965-98 requirements

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